

## Ti Star Internal Fixture System

### Ref. Code And Description

\* (Fixture)


Fixture	Ref. Code	Part Dimension [Unit: mm]	
		Diameter	Length
	FEXIOR3507	3.5	7.0
	FEXIOR3508	3.5	8.5
	FEXIOR3510	3.5	10.0
	FEXIOR3511	3.5	11.5
	FEXIOR3513	3.5	13.0
	FEXIOR3515	3.5	15.0
	FEXIOR4107	4.1	7.0
	FEXIOR4108	4.1	8.5
	FEXIOR4110	4.1	10.0
	FEXIOR4111	4.1	11.5
	FEXIOR4113	4.1	13.0
	FEXIOR4115	4.1	15.0
	FEXIOR4807	4.8	7.0
	FEXIOR4808	4.8	8.5
	FEXIOR4810	4.8	10.0
	FEXIOR4811	4.8	11.5
	FEXIOR4813	4.8	13.0
	FEXIOR4815	4.8	15.0
	FEXIOR5507	5.5	7.0
	FEXIOR5508	5.5	8.5
	FEXIOR5510	5.5	10.0
	FEXIOR5511	5.5	11.5
	FEXIOR5513	5.5	13.0
	FEXIOR5515	5.5	15.0
	FEXIOR6007	6.0	7.0
	FEXIOR6008	6.0	8.5
	FEXIOR6010	6.0	10.0
	FEXIOR6011	6.0	11.5
	FEXIOR6013	6.0	13.0
	FEXIOR6015	6.0	15.0
	FEXIOR7007	7.0	7.0
	FEXIOR7008	7.0	8.5
	FEXIOR7010	7.0	10.0
	FEXIOR7011	7.0	11.5
	FEXIOR7013	7.0	13.0
	FEXIOR7015	7.0	15.0

	FEXIORA3507	3.50	7.0
	FEXIORA3508	3.50	8.5
	FEXIORA3510	3.50	10.0
	FEXIORA3511	3.50	11.5
	FEXIORA3513	3.50	13.0
	FEXIORA3515	3.50	15.0
	FEXIORA4107	4.10	7.0
	FEXIORA4108	4.10	8.5
	FEXIORA4110	4.10	10.0
	FEXIORA4111	4.10	11.5
	FEXIORA4113	4.10	13.0
	FEXIORA4115	4.10	15.0
	FEXIORA4807	4.80	7.0
	FEXIORA4808	4.80	8.5
	FEXIORA4810	4.80	10.0
	FEXIORA4811	4.80	11.5
	FEXIORA4813	4.80	13.0
	FEXIORA4815	4.80	15.0
	FEXIORA5507	5.50	7.0
	FEXIORA5508	5.50	8.5
	FEXIORA5510	5.50	10.0
	FEXIORA5511	5.50	11.5
	FEXIORA5513	5.50	13.0
	FEXIORA5515	5.50	15.0
	FEXIORA6007	6.00	7.0
	FEXIORA6008	6.00	8.5
	FEXIORA6010	6.00	10.0
	FEXIORA6011	6.00	11.5
	FEXIORA6013	6.00	13.0
	FEXIORA6015	6.00	15.0
	FEXIORA7007	7.00	7.0
	FEXIORA7008	7.00	8.5
	FEXIORA7010	7.00	10.0
	FEXIORA7011	7.00	11.5
	FEXIORA7013	7.00	13.0
	FEXIORA7015	7.00	15.0


### Implant Finished Product Packing Contents:

1. With Healing Screw : Ref. Code- IRH503°
2. Without Healing Screw°



\* (Cover Screw) - Material : Ti6Al4V


Cover Screw	Ref. Code	Part Dimension [Unit: mm]	
		Diameter	Post H
	IRC100	5.5	7.0

\* (Closing Screw) - Material : Ti6Al4V



Closing Screw	Ref. Code	Part Dimension [Unit: mm]	
		Diameter	Post H
	ICS200	3.5	6.0

\* (Healing Screw) - Material : Ti6Al4V




Healing Screw	Ref. Code	Part Dimension [Unit: mm]	
		Diameter	Post H
	IRH502	5.5	2.0
	IRH503	5.5	3.0
	IRH504	5.5	4.0
	IRH505	5.5	5.0
	IRH506	5.5	6.0
	IRH507	5.5	7.0
		ATSH503	5.7
ATSH505		5.7	5.3
ATSH507		5.7	7.3
ATSH603		6.7	3.3
ATSH605		6.7	5.3
ATSH607		6.7	7.3


Healing Screw	Ref. Code	Part Dimension [Unit: mm]	
		Diameter	Post H
	ATSH802	8.0	2.3
	ATSH803	8.0	3.3
	ATSH804	8.0	4.3
	ATSH805	8.0	5.3

\* (Abutment Screw) - Material : Ti6Al4V


Abutment Screw	Ref. Code	Part Dimension [Unit: mm]	
		Diameter	Post H
	IAAS600	2.6	5.5
	IRCS200	2.5	4.85
	EXIMS100	2.5	8.6

\* (EZ Post Abutment) - Material : Ti


EZ Post Abutment	Ref. Code	Part Dimension [Unit: mm]		
		Diameter	cuff	Post H
	SEPNOC1	5.5	1.0	7.0
	SEPNOC2	5.5	2.0	7.0
	SEPNOC3	5.5	3.0	7.0
	SEPNOC4	5.5	4.0	7.0
	SEPOCC1	5.5	1.0	7.0
	SEPOCC2	5.5	2.0	7.0
	SEPOCC3	5.5	3.0	7.0
	SEPOCC4	5.5	4.0	7.0
	AEXEPN5515	5.5	1.0	5.5
	AEXEPN5525	5.5	2.0	5.5
	AEXEPN5535	5.5	3.0	5.5
	AEXEPN5545	5.5	4.0	5.5
	AEXEPN6515	6.5	1.0	5.5
	AEXEPN6525	6.5	2.0	5.5
	AEXEPN6535	6.5	3.0	5.5
	AEXEPN6545	6.5	4.0	5.5

EZ Post Abutment	Ref. Code	Part Dimension [Unit: mm]		
		Diameter	cuff	Post H
	ATSSPO5515	5.5	1.5	5.5
	ATSSPO5525	5.5	2.5	5.5
	ATSSPO5535	5.5	3.5	5.5
	ATSSPO5545	5.5	4.5	5.5


\* (Solid Post Abutment) - Material : Ti6Al4V

Solid Post Abutment	Ref. Code	Part Dimension [Unit: mm]		
		Diameter	Cuff	Post H
	ISPA401	4.928	1.0	4.0
	ISPA402	4.928	2.0	4.0
	ISPA403	4.928	3.0	4.0
	ISPA551	4.928	1.0	5.5
	ISPA552	4.928	2.0	5.5
	ISPA553	4.928	3.0	5.5
	ISPA701	4.928	1.0	7.0
	ISPA702	4.928	2.0	7.0
	ISPA703	4.928	3.0	7.0


\* (Solid Abutment) -Material :Ti6Al4V

Solid Abutment	Ref. Code	Part Dimension [Unit: mm]	
		Diameter	Post H
	ISA400	3.5	4.0
	ISA550	3.5	5.5
	ISA700	3.5	7.0



\* (Ball Abutment)- Material :Ti6Al4V; Suitable For Over Denture Case.

Ball Abutment	Ref. Code	Part Dimension [Unit: mm]		
		Diameter	Cuff	Post H
	ATSBA3504	3.5	0.0	4.0
	ATSBA3514	3.5	1.0	4.0
	ATSBA3524	3.5	2.0	4.0
	ATSBA3534	3.5	3.0	4.0
	ATSBA3544	3.5	4.0	4.0


\* (Milling Abutment) -Material :Ti

Milling Abutment	Ref. Code	Part Dimension [Unit: mm]		
		Diameter	cuff	Post H
	ATSP5190	4.9	1.0	9.0
	ATSP5290	4.9	2.0	9.0
	ATSP5390	4.9	3.0	9.0
	ATSP5490	4.9	4.0	9.0


\* (Gold Abutment) -Material :Ti



Gold Abutment	Ref. Code	Part Dimension [Unit: mm]	
		Diameter	Post H
	ATSGA100	5.5	5.65
	ATSGA200	5.5	5.65

\* (Zirconia Abutment) -Material :Ti



Zirconia Abutment	Ref. Code	Part Dimension [Unit: mm]	
		Diameter	Post H
	ATSHA4105	4.1	5.0

\* (Multi Post Abutment & Screw) -Material :Ti6Al4V

Multi Post Abutment & Screw	Ref. Code	Part Dimension [Unit: mm]	
		Diameter	Post H
	EXIM001	5.5	8.7
	EXIM002	5.5	9.7
	EXIM003	5.5	10.7
	EXIM004	5.5	11.7

Multi Post Abutment & Screw	Ref. Code	Part Dimension [Unit: mm]	
		Diameter	Post H
	EXIM651	6.5	8.7
	EXIM652	6.5	9.7
	EXIM653	6.5	10.7
	EXIM654	6.5	11.7
	SMFS002	2.9	11.0

\* (Angled Abutment) -Material :Ti

Angled Abutment	Ref. Code	Part Dimension [Unit: mm]		
		Diameter	Post H	Angle
	IAA3515	3.5	2.5	15°
	IAA3525	3.5	2.5	25°

T-Plus Implant Tech. Co., Ltd.

No.41, Wuquan 6th Rd., Wugu Dist., New Taipei City 248, Taiwan (R.O.C.)

Tel. +886-2-2298-1950 / Fax. +886-2-2298-4353